

10 kDa heat shock protein, mitochondrial

Catalog No: #AP81710



Package Size: #AP81710-1 50ug #AP81710-2 100ug #AP81710-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	10 kDa heat shock protein, mitochondrial
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	Hsp10, 10 kDa chaperonin, Chaperonin 10, CPN10
Accession No.	Uniprot ID:P61603
Uniprot	P61603
GeneID	281833;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: HSPE1

Note: This product is for in vitro research use only