

CAMPATH-1 antigen

Catalog No: #AP81715



Package Size: #AP81715-1 50ug #AP81715-2 100ug #AP81715-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CAMPATH-1 antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Epididymal secretory protein E5,CDW52,CE5,CD52
Accession No.	Uniprot ID:Q28896
Uniprot	Q28896
GeneID	403918;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD52

Note: This product is for in vitro research use only