

CD44 antigen

Catalog No: #AP81718

Package Size: #AP81718-1 50ug #AP81718-2 100ug #AP81718-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD44 antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Extracellular matrix receptor III,ECMR-III,GP90 lymphocyte homing/adhesion receptor,HUTCH-I,Hermes antigen,Hyaluronate receptor,Lymphocyte antigen 24,Ly-24,Phagocytic glycoprotein 1,PGP-1,Phagocytic glycoprotein I,PGP-I,CD44,Ly-24
Accession No.	Uniprot ID:P15379
Uniprot	P15379
GeneID	12505;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cd44

Note: This product is for in vitro research use only