

DNA replication factor Cdt1

Catalog No: #AP81720



Package Size: #AP81720-1 50ug #AP81720-2 100ug #AP81720-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	DNA replication factor Cdt1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Double parked homolog,DUP
Accession No.	Uniprot ID:Q9H211
Uniprot	Q9H211
GeneID	81620;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CDT1

Note: This product is for in vitro research use only