

CUB domain-containing protein 1

Catalog No: #AP81726



Package Size: #AP81726-1 50ug #AP81726-2 100ug #AP81726-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CUB domain-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Membrane glycoprotein gp140, Transmembrane and associated with src kinases, CD318
Accession No.	Uniprot ID:Q5U462
Uniprot	Q5U462
GeneID	109332;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cdcp1

Note: This product is for in vitro research use only