

Cyclin-dependent kinase 12

Catalog No: #AP81727



Package Size: #AP81727-1 50ug #AP81727-2 100ug #AP81727-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin-dependent kinase 12
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Cdc2-related kinase, arginine/serine-rich, CrkRS, Cell division cycle 2-related protein kinase 7, CDC2-related protein kinase 7, Cell division protein kinase 12, hCDK12, CRK7, CRKRS, KIAA0904
Accession No.	Uniprot ID:Q9NYV4
Uniprot	Q9NYV4
GeneID	51755;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CDK12

Note: This product is for in vitro research use only