

Cholinesterase

Catalog No: #AP81730



Package Size: #AP81730-1 50ug #AP81730-2 100ug #AP81730-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cholinesterase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Acylcholine acylhydrolase,Butyrylcholine esterase,Choline esterase II,Pseudocholinesterase,CHE1
Accession No.	Uniprot ID:P06276
Uniprot	P06276
GeneID	590;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:BCHE

Note: This product is for in vitro research use only