

Cyclin-dependent kinase 1

Catalog No: #AP81731



Package Size: #AP81731-1 50ug #AP81731-2 100ug #AP81731-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cyclin-dependent kinase 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CDK1, Cell division control protein 2 homolog, Cell division protein kinase 1, p34 protein kinase, Cdc2, Cdc2a, Cdkn1
Accession No.	Uniprot ID:P39951
Uniprot	P39951
GeneID	54237;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cdk1

Note: This product is for in vitro research use only