

Centromere protein R

Catalog No: #AP81734

Package Size: #AP81734-1 50ug #AP81734-2 100ug #AP81734-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Centromere protein R
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CENP-R,Nuclear receptor-interacting factor 3,Cenpr,Nrif3
Accession No.	Uniprot ID:Q9CQ82
Uniprot	Q9CQ82
GeneID	67733;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Itgb3bp

Note: This product is for in vitro research use only