

C-C motif chemokine 19

Catalog No: #AP81736

Package Size: #AP81736-1 50ug #AP81736-2 100ug #AP81736-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 19
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Epstein-Barr virus-induced molecule 1 ligand chemokine,EBI1 ligand chemokine,Small-inducible cytokine A19,Elc,Scya19
Accession No.	Uniprot ID:O70460
Uniprot	O70460
GeneID	100039053;24047;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ccl19

Note: This product is for in vitro research use only