

C-C motif chemokine 26

Catalog No: #AP81737

Package Size: #AP81737-1 50ug #AP81737-2 100ug #AP81737-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C-C motif chemokine 26
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CC chemokine IMAC, Eotaxin-3, Macrophage inflammatory protein 4-alpha, MIP-4-alpha, Small-inducible cytokine A26, Thymic stroma chemokine-1, TSC-1, SCYA26, UNQ216/PRO242
Accession No.	Uniprot ID:Q9Y258
Uniprot	Q9Y258
GeneID	10344;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CCL26

Note: This product is for in vitro research use only