

# C-C motif chemokine 5

Catalog No: #AP81738



Package Size: #AP81738-1 50ug #AP81738-2 100ug #AP81738-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	C-C motif chemokine 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	SIS-delta, Small-inducible cytokine A5, T-cell-specific protein RANTES, Scya5
Accession No.	Uniprot ID:P50231
Uniprot	P50231
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Ccl5

Note: This product is for in vitro research use only