

CD166 antigen

Catalog No: #AP81739



Package Size: #AP81739-1 50ug #AP81739-2 100ug #AP81739-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD166 antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Activated leukocyte cell adhesion molecule, HB2, KG-CAM, Protein MEMD, SB-10 antigen, CD166
Accession No.	Uniprot ID:O35112
Uniprot	O35112
GeneID	79559;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Alcam

Note: This product is for in vitro research use only