

Myeloid cell surface antigen CD33

Catalog No: #AP81740



Package Size: #AP81740-1 50ug #AP81740-2 100ug #AP81740-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Myeloid cell surface antigen CD33
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Sialic acid-binding Ig-like lectin 3, Siglec-3, CD33, Siglec3
Accession No.	Uniprot ID:Q63994
Uniprot	Q63994
GeneID	12489;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cd33

Note: This product is for in vitro research use only