

# CD276 antigen

Catalog No: #AP81744



Package Size: #AP81744-1 50ug #AP81744-2 100ug #AP81744-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CD276 antigen
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	B7 homolog 3,B7-H3, Costimulatory molecule,CD276,B7h3
Accession No.	Uniprot ID:Q8VE98
Uniprot	Q8VE98
GeneID	102657;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Cd276

Note: This product is for in vitro research use only