

## Chymase

Catalog No: #AP81747



Package Size: #AP81747-1 50ug #AP81747-2 100ug #AP81747-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Chymase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Alpha-chymase,Mast cell protease 3,rMCP-3,Mast cell protease 5,rMCP-5,Mast cell protease III,rMCP-III,Mcpt3
Accession No.	Uniprot ID:P50339
Uniprot	P50339
GeneID	25627;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cma1

Note: This product is for in vitro research use only