

# Cartilage intermediate layer protein 1

Catalog No: #AP81751



Package Size: #AP81751-1 50ug #AP81751-2 100ug #AP81751-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cartilage intermediate layer protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CILP-1, Cartilage intermediate-layer protein, UNQ602/PRO1188
Accession No.	Uniprot ID:O75339
Uniprot	O75339
GeneID	8483;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: CILP

Note: This product is for in vitro research use only