

Cold-inducible RNA-binding protein

Catalog No: #AP81752

Package Size: #AP81752-1 50ug #AP81752-2 100ug #AP81752-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Cold-inducible RNA-binding protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	A18 hnRNP, Glycine-rich RNA-binding protein CIRP, Cirp
Accession No.	Uniprot ID:P60824
Uniprot	P60824
GeneID	12696;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cirbp

Note: This product is for in vitro research use only