

## Centromere protein R

Catalog No: #AP81753



Package Size: #AP81753-1 50ug #AP81753-2 100ug #AP81753-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Centromere protein R
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CENP-R,Beta-3-endonexin,Integrin beta-3-binding protein,Nuclear receptor-interacting factor 3,CENPR,NRIF3
Accession No.	Uniprot ID:Q13352
Uniprot	Q13352
GeneID	23421;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ITGB3BP

Note: This product is for in vitro research use only