

Complement factor I

Catalog No: #AP81754

Package Size: #AP81754-1 50ug #AP81754-2 100ug #AP81754-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Complement factor I
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	C3B/C4B inactivator,If
Accession No.	Uniprot ID:Q61129
Uniprot	Q61129
GeneID	12630;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cfi

Note: This product is for in vitro research use only