

T-lymphocyte activation antigen CD80

Catalog No: #AP81755

Package Size: #AP81755-1 50ug #AP81755-2 100ug #AP81755-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	T-lymphocyte activation antigen CD80
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Activation B7-1 antigen,B7,CD80,B7
Accession No.	Uniprot ID:Q00609
Uniprot	Q00609
GeneID	12519;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name: Cd80

Note: This product is for in vitro research use only