

Eotaxin

Catalog No: #AP81758



Package Size: #AP81758-1 50ug #AP81758-2 100ug #AP81758-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Eotaxin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	C-C motif chemokine 11, Small-inducible cytokine A11
Accession No.	Uniprot ID:Q9TTS6
Uniprot	Q9TTS6
GeneID	404072;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CCL11

Note: This product is for in vitro research use only