

Serine/threonine-protein kinase Chk1

Catalog No: #AP81759



Package Size: #AP81759-1 50ug #AP81759-2 100ug #AP81759-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine/threonine-protein kinase Chk1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	CHK1 checkpoint homolog,Checkpoint kinase-1,CHK1
Accession No.	Uniprot ID:Q8AYC9
Uniprot	Q8AYC9
GeneID	374260;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CHEK1

Note: This product is for in vitro research use only