

DNA replication factor Cdt1

Catalog No: #AP81761



Package Size: #AP81761-1 50ug #AP81761-2 100ug #AP81761-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DNA replication factor Cdt1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Double parked homolog,DUP,Retroviral insertion site 2 protein,Ris2
Accession No.	Uniprot ID:Q8R4E9
Uniprot	Q8R4E9
GeneID	67177;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cdt1

Note: This product is for in vitro research use only