

Calcium and integrin-binding protein 1

Catalog No: #AP81764



Package Size: #AP81764-1 50ug #AP81764-2 100ug #AP81764-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calcium and integrin-binding protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CIB, Calmyrin, DNA-PKcs-interacting protein, Kinase-interacting protein, KIP, Cib, Kip, Prkdcip
Accession No.	Uniprot ID:Q9Z0F4
Uniprot	Q9Z0F4
GeneID	23991;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cib1

Note: This product is for in vitro research use only