

Cerebral dopamine neurotrophic factor

Catalog No: #AP81770



Package Size: #AP81770-1 50ug #AP81770-2 100ug #AP81770-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Cerebral dopamine neurotrophic factor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	ARMET-like protein 1,Conserved dopamine neurotrophic factor,ARMETL1
Accession No.	Uniprot ID:Q49AH0
Uniprot	Q49AH0
GeneID	441549;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CDNF

Note: This product is for in vitro research use only