

# Ceramide glucosyltransferase

Catalog No: #AP81771



Package Size: #AP81771-1 50ug #AP81771-2 100ug #AP81771-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Ceramide glucosyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	GLCT-1, Glucosylceramide synthase, GCS, UDP-glucose ceramide glucosyltransferase, UDP-glucose:N-acylsphingosine D-glucosyltransferase
Accession No.	Uniprot ID:Q16739
Uniprot	Q16739
GeneID	7357;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: UGCG

Note: This product is for in vitro research use only