

Calcium-binding mitochondrial carrier protein Aralar2

Catalog No: #AP81776



Package Size: #AP81776-1 50ug #AP81776-2 100ug #AP81776-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Calcium-binding mitochondrial carrier protein Aralar2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Citrin, Mitochondrial aspartate glutamate carrier 2, Solute carrier family 25 member 13, Aralar2
Accession No.	Uniprot ID:Q9QXX4
Uniprot	Q9QXX4
GeneID	50799;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Slc25a13

Note: This product is for in vitro research use only