

## CCAAT/enhancer-binding protein beta

Catalog No: #AP81777



Package Size: #AP81777-1 50ug #AP81777-2 100ug #AP81777-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CCAAT/enhancer-binding protein beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	C/EBP beta,C/EBP-related protein 2,Interleukin-6-dependent-binding protein,IL-6DBP,Liver-enriched inhibitory protein,LIP,Liver-enriched transcriptional activator,LAP,Silencer factor B,SF-B,Crp2,Nf- $\kappa$ B,Sfb
Accession No.	Uniprot ID:P21272
Uniprot	P21272
GeneID	24253;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cebpb

Note: This product is for in vitro research use only