

Chitinase-3-like protein 1

Catalog No: #AP81779



Package Size: #AP81779-1 50ug #AP81779-2 100ug #AP81779-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chitinase-3-like protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	38 kDa heparin-binding glycoprotein,Signal-processing protein,gp38k
Accession No.	Uniprot ID:Q29411
Uniprot	Q29411
GeneID	396865;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CHI3L1

Note: This product is for in vitro research use only