

Carcinoembryonic antigen-related cell adhesion molecule 1

Catalog No: #AP81785

Package Size: #AP81785-1 50ug #AP81785-2 100ug #AP81785-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Carcinoembryonic antigen-related cell adhesion molecule 1 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Hu |
| Other Names | Biliary glycoprotein 1,BGP-1,CD66a,BGP,BGP1 |
| Accession No. | Uniprot ID:P13688 |
| Uniprot | P13688 |
| GeneID | 634; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CEACAM1

Note: This product is for in vitro research use only