

Choline kinase alpha

Catalog No: #AP81789

Package Size: #AP81789-1 50ug #AP81789-2 100ug #AP81789-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Choline kinase alpha
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	CK,CHETK-alpha,Ethanolamine kinase,EK,Chk,Ckr
Accession No.	Uniprot ID:Q01134
Uniprot	Q01134
GeneID	29194;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Chka

Note: This product is for in vitro research use only