

## CMRF35-like molecule 7

Catalog No: #AP81791

Package Size: #AP81791-1 50ug #AP81791-2 100ug #AP81791-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |  |
|--------------------|--|
| Product Name       | CMRF35-like molecule 7   |
| Brief Description  | Recombinant Protein  |
| Host Species       | E.coli   |
| Purification       | Greater than 90% by SDS-PAGE   |
| Species Reactivity | Hu   |
| Other Names        | CLM-7,CD300 antigen-like family member B,CMRF35-A2,Immune receptor expressed on myeloid cells 3,IREM-3,Leukocyte mono-Ig-like receptor 5,Triggering receptor expressed on myeloid cells 5,TREM-5,CD300b,CD300B,CLM7,CMRF35A2,IREM3,LMIR5,TREM5,UNQ2530/PRO6029 |
| Accession No.      | Uniprot ID:A8K4G0  |
| Uniprot            | A8K4G0   |
| GeneID             | 124599;  |
| Tag Info           | His  |
| Storage            | Store at -20°C. (Avoid repeated freezing and thawing.)   |

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CD300LB

Note: This product is for in vitro research use only