

C-C motif chemokine 13

Catalog No: #AP81792



Package Size: #AP81792-1 50ug #AP81792-2 100ug #AP81792-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | C-C motif chemokine 13 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Dog |
| Other Names | Monocyte chemoattractant protein 4, Monocyte chemotactic protein 4, MCP-4, Small-inducible cytokine A13 |
| Accession No. | Uniprot ID:Q68Y88 |
| Uniprot | Q68Y88 |
| GeneID | 445451; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CCL13

Note: This product is for in vitro research use only