

## C-C motif chemokine 21

Catalog No: #AP81793



Package Size: #AP81793-1 50ug #AP81793-2 100ug #AP81793-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	C-C motif chemokine 21
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	6Ckine,Beta-chemokine exodus-2,Secondary lymphoid-tissue chemokine,SLC,Small-inducible cytokine A21,SCYA21,UNQ784/PRO1600
Accession No.	Uniprot ID:O00585
Uniprot	O00585
GeneID	6366;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CCL21

Note: This product is for in vitro research use only