

## Cytosolic non-specific dipeptidase

Catalog No: #AP81795



Package Size: #AP81795-1 50ug #AP81795-2 100ug #AP81795-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cytosolic non-specific dipeptidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	CNDP dipeptidase 2, Glutamate carboxypeptidase-like protein 1, Cn2
Accession No.	Uniprot ID:Q9D1A2
Uniprot	Q9D1A2
GeneID	66054;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Cndp2

Note: This product is for in vitro research use only