

## Chymase

Catalog No: #AP81797



Package Size: #AP81797-1 50ug #AP81797-2 100ug #AP81797-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Chymase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Alpha-chymase,Mast cell chymase 1,Mast cell protease 5,mMCP-5,Mast cell protease I,Mcpt5
Accession No.	Uniprot ID:P21844
Uniprot	P21844
GeneID	17228;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Cma1

Note: This product is for in vitro research use only