

Carcinoembryonic antigen-related cell adhesion molecule 5

Catalog No: #AP81798

Package Size: #AP81798-1 50ug #AP81798-2 100ug #AP81798-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Carcinoembryonic antigen-related cell adhesion molecule 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Pregnancy-specific glycoprotein 30,Psg30
Accession No.	Uniprot ID:Q3UUK2
Uniprot	Q3UUK2
GeneID	73250;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ceacam5

Note: This product is for in vitro research use only