

Chondroitin sulfate proteoglycan 4

Catalog No: #AP81800



Package Size: #AP81800-1 50ug #AP81800-2 100ug #AP81800-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Chondroitin sulfate proteoglycan 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Chondroitin sulfate proteoglycan NG2, HSN tumor-specific antigen, Ng2
Accession No.	Uniprot ID:Q00657
Uniprot	Q00657
GeneID	81651;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Cspg4

Note: This product is for in vitro research use only