

Chymase

Catalog No: #AP81801



Package Size: #AP81801-1 50ug #AP81801-2 100ug #AP81801-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Chymase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Dog
Other Names	Alpha-chymase,Mast cell protease I
Accession No.	Uniprot ID:P21842
Uniprot	P21842
GeneID	490628;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CMA1

Note: This product is for in vitro research use only