

Claudin-3

Catalog No: #AP81804



Package Size: #AP81804-1 50ug #AP81804-2 100ug #AP81804-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Claudin-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Clostridium perfringens enterotoxin receptor 2,CPE-R 2,Cpetr2
Accession No.	Uniprot ID:Q9Z0G9
Uniprot	Q9Z0G9
GeneID	12739;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Cldn3

Note: This product is for in vitro research use only