

T-cell-specific surface glycoprotein CD28 homolog

Catalog No: #AP81805



Package Size: #AP81805-1 50ug #AP81805-2 100ug #AP81805-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	T-cell-specific surface glycoprotein CD28 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	CHT28
Accession No.	Uniprot ID:P31043
Uniprot	P31043
GeneID	396249;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD28

Note: This product is for in vitro research use only