

T-lymphocyte activation antigen CD86

Catalog No: #AP81806



Package Size: #AP81806-1 50ug #AP81806-2 100ug #AP81806-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	T-lymphocyte activation antigen CD86
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Activation B7-2 antigen,B70,BU63,CTLA-4 counter-receptor B7.2,FUN-1,CD86,CD28LG2
Accession No.	Uniprot ID:P42081
Uniprot	P42081
GeneID	942;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:CD86

Note: This product is for in vitro research use only