

Collagen alpha-1(VII) chain

Catalog No: #AP81808



Package Size: #AP81808-1 50ug #AP81808-2 100ug #AP81808-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Collagen alpha-1(VII) chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Long-chain collagen,LC collagen
Accession No.	Uniprot ID:Q02388
Uniprot	Q02388
GeneID	1294;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:COL7A1

Note: This product is for in vitro research use only