

Complement component C8 beta chain

Catalog No: #AP81809



Package Size: #AP81809-1 50ug #AP81809-2 100ug #AP81809-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Complement component C8 beta chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Complement component 8 subunit beta
Accession No.	Uniprot ID:P55314
Uniprot	P55314
GeneID	313421;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: C8b

Note: This product is for in vitro research use only