

STIP1 homology and U box-containing protein 1

Catalog No: #AP81811



Package Size: #AP81811-1 50ug #AP81811-2 100ug #AP81811-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	STIP1 homology and U box-containing protein 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Carboxy terminus of Hsp70-interacting protein,E3 ubiquitin-protein ligase CHIP,RING-type E3 ubiquitin transferase CHIP,Chip
Accession No.	Uniprot ID:Q9WUD1
Uniprot	Q9WUD1
GeneID	56424;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Stub1

Note: This product is for in vitro research use only