

Choline O-acetyltransferase

Catalog No: #AP81812



Package Size: #AP81812-1 50ug #AP81812-2 100ug #AP81812-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Choline O-acetyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Pig
Other Names	CHOACTase
Accession No.	Uniprot ID:P13222
Uniprot	P13222
GeneID	396896;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: CHAT

Note: This product is for in vitro research use only