

## Serine/threonine-protein kinase Chk1

Catalog No: #AP81816



Package Size: #AP81816-1 50ug #AP81816-2 100ug #AP81816-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Serine/threonine-protein kinase Chk1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CHK1 checkpoint homolog,Cell cycle checkpoint kinase,Checkpoint kinase-1,CHK1
Accession No.	Uniprot ID:O14757
Uniprot	O14757
GeneID	1111;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:CHEK1

Note: This product is for in vitro research use only