

Clusterin

Catalog No: #AP81821



Package Size: #AP81821-1 50ug #AP81821-2 100ug #AP81821-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Clusterin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Apolipoprotein J,Apo-J,Clustrin,Sulfated glycoprotein 2,SGP-2,Apoj,Msgp-2
Accession No.	Uniprot ID:Q06890
Uniprot	Q06890
GeneID	12759;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Clu

Note: This product is for in vitro research use only