

Complement component C8 alpha chain

Catalog No: #AP81824



Package Size: #AP81824-1 50ug #AP81824-2 100ug #AP81824-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Complement component C8 alpha chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Complement component 8 subunit alpha
Accession No.	Uniprot ID:Q8K182
Uniprot	Q8K182
GeneID	230558;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: C8a

Note: This product is for in vitro research use only